



## The PHENIX CENTRAL MAGNET OPERATION INSTRUCTIONS

procedure name

**PHENIX Procedure No. PP-2.5.1.3-01**

**Revision: C**

**Date: 1/25/2013**

### **Hand Processed Changes**

**HPC No.**

**Date**



**Page Nos.**

**Initials**

- Typo: Under 3.0 Prerequisites, PHENIX Procedure No. PP-2.5.1.3-01, "Procedure for Preparation of Central Magnet for Operation" should be written as PHENIX Procedure No. PP-2.5.1.3-02, "Procedure for Preparation of Central Magnet for Routine Operation"*

_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____

### **Approvals**

	
PHENIX S E & I	Cognizant Scientist/Engineer
Date	Date
	/Activity Manager

	1-30-13
PHENIX QA/Safety	Date

*The only official copy of this file is the one online. Before using a printed copy, verify that it is the most current version by checking the document effective date on this web site.*



# The PHENIX CENTRAL MAGNET OPERATION INSTRUCTIONS

---

procedure name

**PHENIX Procedure No. PP-2.5.1.3-01**

**Revision: C**

**Date: 1/25/2013**

## **Hand Processed Changes**

<b><u>HPC No.</u></b>	<b><u>Date</u></b>	<b><u>Page Nos.</u></b>	<b><u>Initials</u></b>
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____

## **Approvals**

_____ PHENIX S E & I	_____ Date	_____ Cognizant Scientist/Engineer /Activity Manager	_____ Date
_____ PHENIX QA/Safety	_____ Date	_____	

## REVISION CONTROL SHEET

LETTER	DESCRIPTION	DATE	WRITTEN BY	APPROVED BY	Current Oversight
A	First Issue	03/02/2000	n/a	J. Haggerty, W. Lenz & (unintelligible)	n/a
B	Reviewed and found to be OK as is	11/13/2009	D. Lynch	P. Giannotti, D. Lynch, R. Pisani	P. Giannotti
C	Reviewed and found to be OK as is	1/25/2013	D. Lynch	P. Giannotti, D. Lynch, R. Pisani	P. Giannotti

**Procedure for Turning on the PHENIX Central Magnet**

**1.0 Purpose and Scope**

The purpose of this procedure is to turn on and ramp the PHENIX Central Magnet (CM) to the desired field.

**2.0 Responsibilities**

The PHENIX shift leader is responsible for conducting the procedure.

**3.0 Prerequisites**

The PHENIX Run Coordinator shall have completed the requirements of PHENIX Procedure No. PP-2.5.1.3-01, "Procedure for Preparation of Central Magnet for Operation."

**4.0 Precautions**

None

**5.0 Procedure**

5.1 The PHENIX shift leader shall call the main control room to ask for approval to excite or change the settings of the PHENIX magnet. The power level should be given, namely 700 kW for full field or 175 kW for half field.

5.2 The PHENIX shift leader shall then use the PHENIX magnet control page to set the ramp speed followed by the current values for the magnet power supply. Copies of the control page giving the settings for full and half field have been given to all PHENIX Magnet Operators.

5.3 The PHENIX shift leader shall then excite the coils.

**6.0 Documentation**

None

**7.0    References**

None

**8.0    Attachments**

None