

Running the Data Acquisition

Run 03

B. Cole, S. Adler, J. Haggerty, J. M.
Burward-Hoy, S. Batsouli

Master Clock


1. Make sure you are on the two-panel display, labeled “phoncsb”.
2. Click on the Gnome paw on the lower left hand corner to open the ONCS Window. This is the menu of all the functions we will use for the data acquisition.
3. Choose the “Master Clock” gui
4. Select the type of clock depending on what we have to do:

BEAM is COGGED – SL will tell you to run the DAQ.

- Select the “Blue Clock” tab. The window should then say “Clock Source is Blue Clock.”

WHEN SL SAYS THAT MCR WILL DUMP THE BEAM



- Run “Master Clock” gui and select the “Internal Clock” tab. The window should then say “Clock Source is Internal Clock”.



Global Level 1 (GL1) Trigger



Go to ONCS Window Menu

- Run “Trigger” gui.
- On GL1 Tab, select “Run Config. Script”. A smaller window menu will pop up asking you to select the type of trigger. This depends on the clock source as follows:
 - If **Blue Clock**: Choose “BigRun3”, then “Load”
 - OR
 - If **Internal Clock**: Choose the configuration script depending on subsystem requests- ask SL, then “Load”
- Select BBC Tab
 - Select “Run Config. Script”.
 - For either blue or internal clocks, select “PP” (for BBLL1).





Parallel Clock Fix on the FEED Gui

1. On the ONCS Menu, select “feed”.
Feed stand for “Front End Electronics Defibrillator.” This initializes the FEMs so that the phase-lock loops are in phase in the DLINKs.
2. We’ve changed the clock, so select “Parallel Clock Fix.”
3. After the green progress bar finishes, a green status report will appear below the “Parallel Clock Fix” button.
4. If there is a problem with this initialization for a given subsystem, its status indicator will be red. If this is the case, discuss this with the SL if this is the normal state for the detector. If not, push the button for the given subsystem to initialize again. If it continues to fail check with SA2 about LV status and consult SL



Monitoring the Crates on PHONCS C


- Go to the computer labeled “PHONCS C”.
- Select “PPC Windows” on the CNCS menu.
- The green windows are the GTM crates.
- The pink windows are GL1 and BBC LL1 crates.
- The blue windows are the iocondev crates for the granules.
- These windows are useful for monitoring during the download step in run control (error messages or “hanging” prompts will occur if there are problems during run control setup).



Monitoring the Crates on PHONCS C_{cont}

- If “hanging” prompts occur (i.e. memory allocation error) then reboot the relevant iocondev by pressing Ctrl-x
- If rebooting does not help or if the error “event should be 1 instead of 0” appears contact the daq expert (the power might need to be recycled)

Starting the Servers





Event Builder: choose "Enable" in "Data Logging Control." Then choose "Initialize EvB".

Partition Server: select "on" in "Use EvB".

Any errors? The "ok" status button may turn red and error messages will appear in the white scroll bar.


Any errors in the Event Builder? The seb/atp buttons may turn red during run (see for example, ATP Read Errors) and messages will display on the white scroll bar.

If after initial evb initialization the buttons of certain seb/atp are red try clicking on them to reinitialize them individually. If that doesn't work contact the daq expert



Global Level 1 Run Control



- On the ONCS Menu, select “Run Control”. Type your name in the prompt for the name. The name will appear in the label for run control window.
- Select “GL1” from the pop-up window. A run control window will now appear, with the name “Run Control for GL1 Owned by ...”. All the run control commands for this window will appear in the upper left-hand corner of the window.
- In this menu, under “Mode”, select “Evb Enable”.
- Select the “Download” button. This will take approximately one minute.



- If everything is set properly, you should start to see the EVB Status counters increasing. This is where the mean event rate is displayed.
- The trigger status at the bottom should say “Enabled” with raw, live, etc. trigger rates displayed as non-zero entries.

Big Partition Run Control

- On the ONCS Menu, select “Run Control”. Type your name in the prompt for the name. The name will appear in the label for run control window.
- Select “Big” from the pop-up window. A run control window will now appear, with the name “Run Control for BigRun3 Owned by ...”. All the run control commands will appear in the upper left-hand corner of the window.
- In this menu, under “Mode”, select “Evb Enable”. Then select “Select Trigger Manually” according to the SL instructions on which triggers to run.
- Select the “Download” button. This will take approximately one minute. Monitor the crate windows on the PHONCS C terminal during this procedure for any error messages. Now press “Start”. A window will appear. Choose “Physics”.




The trigger status at the bottom should say “Enabled” with raw, live, etc. trigger rates displayed as non-zero entries.

If everything is set properly, you should start to see the EVB Status counters increasing. This is where the mean event rate is displayed.

Big Partition Run Control

- The Trigger settings, scalar values, Scalar monitor should be checked. To do this click at **Config. readout** and **Trigger Settings** from rc (check standing orders with SL)
- Save the settings in the log book. The command for capturing a window is: **“gimp”**, then choose **“acquire”**-> **“screen shot ->”single screen”**
Right-click on the window you want to capture and save.



Available Triggers for Partition Big	
Clock	1,000,000
BBCLL1>=1	15
ZDCNS	16
ERT_2x2	0
ERT_2x2&BBCLL1	4
ERT_Gamma1	0
ERT_Gamma1&BBCLL1	0
ERT_Gamma2	0
ERT_Gamma2&BBCLL1	0
ERT_Electron	0
ERT_Electron&BBCLL1	0
ERT_Gamma1&ERT_Gamma2	0
ERT_Gamma1 ERT_Gamma2	0
ERT_Phi&BBCLL1	4
ERT_Gamma3	0
MUIDS_1D	10,000
MUIDS_1D&BBCLL1	10
MUIDS_1D1S	500
MUIDS_1D1S*BBCLL1	1
MUIDN_1D	10,000
MUIDN_1D&BBCLL1	0
MUIDN_1D1S	100
MUIDN_1D1S*BBCLL1	0
NTCS	0
NTCN	0
NTCNSnarrow	0
NTCNSwide	100
PPG(Pedestal)	0
PPG(Test Pulse)	0
PPG(Laser)	0
BBCLL1>=1(noVertexCut)	200

Set Triggers EXIT

Selected triggers are ENABLED.

ATP Read Errors

- Monitor the “ATP Status”-“Read Err” column. This column should all read “0”.
- If there are numbers, indicating errors, then we must stop the run.
- It may be that the FEED step needs to be repeated for a given subsystem. All subsystems should receive the same number of events, so scan through the column “#Events” and look for that subsystem that has a different event rate from the others. In the example case below, the TEC.E has a different event number. So, the FEED was run for the TEC alone after stopping the Big Partition Run Control.

Granule Names	LT	Run	Busy	OK	DCM Status	Name	Events	LT Size	Data Rate	Buff Usage	Busy	OK	Name	#Events	#LT2Accept	#LT2Gr	ATP Status	Assn Rate	Ave Data Rate	ATP OK	ET OK	EBC Status	EBC.0
BB	657	Run	0	0	OK	SERBB.0	657	552 KB	0.026 MB/s	0.019	0	0	ATP.0	33	0	9	0.678/s	0.075 MB/s	0.075 MB/s	0.075 MB/s	0.075 MB/s	EBC OK	EBC.0
ZDC	657	Run	0	0	OK	SERZDC.0	657	764 KB	0.008 MB/s	0.077	0	0	ATP.1	27	0	16	0.435/s	0.048 MB/s	0.048 MB/s	0.048 MB/s	0.048 MB/s	#Assigned	657
DC.W	657	Run	0	0	OK	SERDC.W.0	657	421 KB	0.050 MB/s	0.077	0	0	ATP.2	27	0	15	0.427/s	0.047 MB/s	0.047 MB/s	0.047 MB/s	0.047 MB/s	#Completed	617
PC.W	657	Run	0	0	OK	SERPC.W.0	657	379 KB	0.040 MB/s	0.077	0	0	ATP.3	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	Avg Event Rate	0.800/s
TEC.E	657	Run	0	0	OK	SERTEC.E.0	657	42 KB	0.049 MB/s	0.077	0	0	ATP.8	28	0	15	0.443/s	0.049 MB/s	0.049 MB/s	0.049 MB/s	0.049 MB/s	Avg Assem Lat	16.424 s
TEC.E	657	Run	0	0	OK	SERTEC.E.1	657	45 KB	0.047 MB/s	0.080	0	0	ATP.9	28	0	15	0.443/s	0.049 MB/s	0.049 MB/s	0.049 MB/s	0.049 MB/s	Avg ATP Load	0.000
TEC.E	657	Run	0	0	OK	SERTEC.E.2	657	64 KB	0.066 MB/s	0.039	0	0	ATP.A	26	0	15	0.412/s	0.046 MB/s	0.046 MB/s	0.046 MB/s	0.046 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.3	657	329 KB	0.013 MB/s	0.087	0	0	ATP.B	27	0	15	0.427/s	0.046 MB/s	0.046 MB/s	0.046 MB/s	0.046 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.4	657	127 KB	0.131 MB/s	0.076	0	0	ATP.C	27	0	15	0.427/s	0.048 MB/s	0.048 MB/s	0.048 MB/s	0.048 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.5	657	3.8 KB	0.039 MB/s	0.076	0	0	ATP.D	30	0	18	0.780/s	0.087 MB/s	0.087 MB/s	0.087 MB/s	0.087 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.6	657	3.5 KB	0.037 MB/s	0.077	0	0	ATP.E	29	0	15	0.843/s	0.094 MB/s	0.094 MB/s	0.094 MB/s	0.094 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.7	657	3.7 KB	0.039 MB/s	0.037	0	0	ATP.F	30	0	18	0.780/s	0.087 MB/s	0.087 MB/s	0.087 MB/s	0.087 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.8	657	3.5 KB	0.037 MB/s	0.077	0	0	ATP.G	29	0	17	0.790/s	0.089 MB/s	0.089 MB/s	0.089 MB/s	0.089 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.9	657	4.0 KB	0.042 MB/s	0.023	0	0	ATP.H	28	0	15	0.443/s	0.049 MB/s	0.049 MB/s	0.049 MB/s	0.049 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.10	657	2.1 KB	0.024 MB/s	0.077	0	0	ATP.I	30	0	18	0.780/s	0.086 MB/s	0.086 MB/s	0.086 MB/s	0.086 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.11	657	4.1 KB	0.043 MB/s	0.154	0	0	ATP.J	30	0	18	0.780/s	0.087 MB/s	0.087 MB/s	0.087 MB/s	0.087 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.12	657	4.7 KB	0.042 MB/s	0.154	0	0	ATP.K	25	0	15	0.390/s	0.045 MB/s	0.045 MB/s	0.045 MB/s	0.045 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.13	657	4.0 KB	0.050 MB/s	0.152	0	0	ATP.L	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.14	657	0.126 KB	0.009 MB/s	0.019	0	0	ATP.M	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.15	657	0.112 KB	0.010 MB/s	0.077	0	0	ATP.N	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.16	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.O	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.17	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.P	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.18	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.Q	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.19	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.R	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.20	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.S	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.21	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.T	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.22	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.U	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.23	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.V	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.24	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.W	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.25	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.X	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.26	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.Y	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.27	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.Z	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.28	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AA	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.29	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AB	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.30	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AC	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.31	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AD	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.32	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AE	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.33	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AF	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.34	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AG	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.35	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AH	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.36	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AI	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.37	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AJ	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.38	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AK	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.39	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AL	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.40	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AM	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.41	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AN	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.42	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AO	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.43	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AP	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.44	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AQ	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.45	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AR	30	0	15	0.475/s	0.053 MB/s	0.053 MB/s	0.053 MB/s	0.053 MB/s		
TEC.E	657	Run	0	0	OK	SERTEC.E.46	657	0.13 KB	0.010 MB/s	0.076	0	0	ATP.AS	30	0	15							

DAQ Monitor for the Run

- Look for GL1 and FEM CLK errors shown as red bars in the picture displayed.
- Make sure the Offline Shift Person monitors the DAQ Monitor for each run.
- If these occur for a given subsystem, and if the SL agrees, stop a run and select the button for the subsystem on the feed gui.
- The message “STOP THE DAQ” at the DAQ monitor means... you should stop the run

