
## Hadron Blind Detector implementation during PHENIX Run-10






On behalf of the PHENIX Collaboration





#### Gas Electron Multiplier



- Invented by F. Sauli at CERN (NIM A386 (1997) 531-534)
- Insulating film (Kapton) sandwiched between layers of metal (Cu, Au)
- HV creates strong field so that an avalanche can occur inside the holes

Gas Gains ~ 10-20 / GEM foil

~ 10<sup>3</sup>-10<sup>4</sup> or higher in triple-GEM configuration




70 um

50 µm

300-500 V


60-100 KV/cm



 Triple GEM stack with wire mesh and Readout Pads








- Triple GEM stack with wire mesh and Readout Pads
- Ionization from charged particles avalanche and charge is collected on Pads






 Apply reverse-bias voltage on mesh-GEM → primary ionization drift away from GEM








- Apply reverse-bias voltage on mesh-GEM → primary ionization drift away from GEM
- Deposit photocathode (CsI) on top GEM → UV photons produce photoelectrons from the CsI photocathode
  - Photoelectrons avalanche in the holes, charge collected by Readout Pads.
  - Triple GEM stack yields a gain of a few x 10<sup>3</sup>





#### The Degree of Hadron Blindness




 At slightly negative E<sub>d</sub>, photoelectron detection efficiency is preserved while charge collection is largely suppressed.







# Proximity Focused Windowless Cherenkov Detector (HBD)





#### HBD in PHENIX






Backgrounds due to low momentum electrons from Dailtz pairs & conversions can be reduced with HBD detection of electrons





## **HBD** History

- Fall 2006 installed in PHENIX for Run-7 (2006-07).
  - HBD's GEM foils were damaged due to severe HV problems
    - Minor GEM sparks induced damaging mesh to GEM sparks.
    - A spark would induce more sparks in other modules due to a copious production of photons from original spark.
- Rebuilt & installed for Run-9 (Feb. July' 09)
  - Built HBD-East using time consuming "test, test, and install method".
  - Built HBD-West using M. Durham's Rapid Assembly Method.







#### HBD History ...

- HBD-East Problems during Run-9:
  - Mesh-GEM short in 1 module (early in run)
    - Disabled module ES1.
  - Another module had deteriorating performance.
  - Sparking resulting in trips cause subpar data collection.
- HBD-West had no major problems during Run-9
- Decision to rebuild HBD-East
  - Known difference in assembly method of East/West arms.
  - Deterioration in performance continuing into Run-10 a possibility.
  - Time scale to rebuild was sufficient.







## Rebuilding of HBD-East

- M. Durham's Rapid Assembly Method.
  - Individual GEMs: clean, condition w/HV, test.
  - Assemble bottom & middle GEMs in an easy to access, clean environment (adjacent to laminar flow table).
  - Install top GEMs (Csl coated) in Nitrogen filled glove box.

Seal HBD vessel, then test GEM stacks.









#### Problem with HBD module EN2


- Testing of HBD-East showed that GEM stack
  EN2 could not hold the nominal operating voltage
  - Each individual GEM able to maintain a voltage exceeding nominal operating voltage.
  - Possibly due to small dent in a GEM.
  - Workaround solution adopted.
    - Revert voltage divider to older design, which maintains a lower voltage
    - Eventually shorts developed leading to trips... EN2 was disabled.
- Other attempts at fixing EN2 have failed
  - Currently EN2 is disabled
  - 5% acceptance loss





## Creating Au+Au Events from p+p events

- HBD is very efficient at determining tracks from p+p events.
- Instead of using Monte Carlo, use real data by accumulating many p+p events to emulate one Au+Au event.
  - p+p tracks for each individual events are well determined.
  - HBD reconstruction will be less efficient in determining Au+Au tracks due to overlapping of tracks.
  - Obtain electron identification in Au+Au collisions by comparing "Accumulated" Au+Au events to known p+p reconstructions.



#### Conclusions

- Rapid Assembly Method proved successful.
- HBD-East rebuild improved performance
  - Module EN2 disabled, 5% loss of acceptance.
  - Remaining HBD-East modules working as expected.
- p+p event accumulator will allow determination of the HBD's electron identification efficiency.
- Due to the HBD, PHENIX has the added benefit of suppressing Dalitz pairs and photon conversions.
- For more on HBD, stick around for Sky Rolnick's Chiral Symmetry Restoration using HBD following this talk.





#### **Extras**





#### **HBD Detector Parameters**

Acceptance  $|\eta| \le 0.45, \, \Delta \phi = 135^{\circ}$ 

GEM size  $(\phi,z)$  23 x 27 cm<sup>2</sup>

Segmentation 26 strips (0.80 x 27 cm)

2 strips (0.65 x 27 cm)

Number of detector modules per arm 10

Frame 5 mm wide, 0.3mm cross

Hexagonal pad size a = 15.6 mm

Number of pads per arm 960

Dead area within central arm acceptance 6%

Radiation length within central arm acceptance box: 0.92%, gas: 0.54%,

preamps+sockets: 0.66%


Total: 2.12%

Weight per arm (including accessories) <10 kg






## HBD Assembly Area (Clean Tent)








#### Reduced Data Format for Run-10



Old 12 sample format Data volume ~ 20-24 KB/ev




New 3 sample format Data volume ~ 7-8 KB/ev





#### Pedestals and Noise








## Gain Determination Using Scintillation in CF<sub>4</sub>

For low multiplicity events, scintillation produces essentially single photoelectrons on each pad

Fit single p.e. distribution to an exponential

Gain ~ 1/slope

Allows pad by pad calibration





Used as a gain monitor throughout the run







## Setting Gains Using Scintillation





Select on most peripheral events where <Npe> << 1











#### West Gains

HBD: Run 304646, Time: Sun Feb 7 06:22:30 2010 BBC Cut 20











#### **East Gains**

HBD: Run 304646, Time: Sun Feb 7 06:22:30 2010 BBC Cut 20











Setting Reverse Bias Voltages



Need to set reverse bias voltage to optimize pe collection efficiency and minimize hadron response

Requires setting voltage between mesh and top GEM to ~ 5V

Applied voltage is ~ 4KV, so this requires high precision (~ 0.1%) and good stability in the HV PS





S.Rolnick 25

