

AEROGEL CHERENKOV COUNTER

PHENIX Focus




Satoshi Takagi

For the High- p_T Upgrade Team


Contents

- Concept
- What is PHENIX Aerogel Cherenkov Counter?
- Mechanical design & Electronics
- Simulation Activity
- Performance of single cell
- Results & Future Plan
- Summary

Currently, PID is ...



Physics Motivations



✓ From basic strategy, it is natural extension for PHENIX to extend its PID in higher p_T region.

✓ Strong motivation given;

- Jet Quenching !?
- Large suppression of pions at high- p_T , while protons show binary scaling!?

- Meson/baryon puzzle?

- Need to extend PID >

5 GeV/c !!


Concept

➤ PID in high p_T region

- Cherenkov Radiation

➤ Cherenkov Radiator

- Low refractive index
- Best index with RICH(CO₂) is $n \sim 1.01$.



Requirements

- | | |
|---|----------------------|
| - Refractive index : $n \sim 1.01$ | - Momentum threshold |
| - Light yield : >10 p.e. | - Resolving power |
| - Uniformity of the light yield : Needed. | - Easy handling |
| - Occupancy in Au+Au collisions : $<10\%$ | - S/N |

Installation Purpose

➤ To enhance the PID capability of PHENIX !!

Momentum [GeV/c]	0.5	1	2	2.5	3	3.7	4	4.2	5	5.5	6	7	~10
π		Yellow	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Red	Red	Red
K		Yellow	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue
p		Yellow	Yellow	Yellow	Yellow	Yellow	Yellow	Blue	Blue	Blue	Blue	Blue	Blue



TOF



RICH



Aerogel (+ TOF or RICH)

Aerogel : ($n=1.011$)

TOF : 100 ps time resolution


RICH : CO_2 , ($n = 1.00041$)

What is Aerogel Counter ??

(I) Outline

➤ Cherenkov Counter (non-ring-imaging type)


- Cherenkov radiator is Silica Aerogel. (MATSUSHITA, SP-12M)
- Photon is detected by 2 PMTs. (HAMAMATSU, R6233)
- All inner surface is covered with DRP Reflector. (Goretex)
- Integration cube for uniformity of light yield. (Air)



What is Aerogel Counter ??

(II) Aerogel Panel

~ Side view ~





~ Top view ~



- ✓ *Direction of each other cell is reversed*
- *Aerogel locates at the same distance from vertex.*

What is Aerogel Counter ??

(III) Where



- 4.5m from vertex
- 4.0m along z direction
- 15 deg. in phi

RED: Aerogel

YELLOW: Integration sphere

GREEN: PMT

(I) Silica Aerogel

➤ Characteristic

- Refractive index $\sim 1.0114 \pm 0.0008$
 - Silica aerogel with **lowest refractive index commercially available !!**
- Density $\sim 40 \text{ mg/cm}^3$
- Transparent for 10mm thickness
 - 64% @ 400nm, 88% @ 550nm
- Hydrophobic
- Long term stability (KEK-Belle)
- Very fragile





**Many many
Aerogel
tiles!!**


Y. Miake

**He is
godfather of aerogel !!**

(II) PMT



PMT R6233-01HA (HAMAMATSU)





PMT R6233-01HA (HAMAMATSU)

- 3-inch diameter **large !!**
- Gain : $> 10^7$ at -1500 V **high gain !!**
- Q.E. : 30 % **high Q.E. !!**
- Dark Current : 2nA **low noise !!**

(III) HV Divider for PMT

- High Gain compared with HAMAMATSU standard
 - Voltage distribution ratio was modified from standard.
- Low Power consumption
 - We do not need cooling.
- Hand made
 - For Thinner material & less space



Mechanical Design

(IV) Mu-metal Shield

➤ Magnetic field (W1 sector)



- B-field at the place where the aerogel counters is installed : **2~8 Gauss**



➤ Mu-metal Shield


- Thickness & Size of the mu-metal shield has been optimized.

- 0.5mm thick, 80mm long is enough !!



Mechanical Design

(V) Reflector Selection




Reflector	Np.e.	Attenuation Length [cm]
Goretex	11.78	6.5 ± 0.3
Lumirror	10.2	6.9 ± 0.4
Tetratex	9.8	6.3 ± 0.4
Tyvek (Double)	9.1	5.8 ± 0.2
Aluminized Mylar	5.3	4.4 ± 0.2
BlackPaper	2.0	2.6 ± 0.1


Goretex is the best!!

Mechanical Design

(VI) Reflector Property



5000x SEM of
DRP® Reflector




- A 3.0 mm DRP® Reflector
- B 0.5 mm DRP® Reflector**
- C 0.25 mm DRP® Reflector
- D Granular PTFE
- E Barium Sulfate
- F Microporous Polyester
- G Powder Coating

Thickness	0.5 mm
Thermal conductivity	0.04 W/m/K
Water resistance	Highly hydrophobic
UV Exposure	UV-resistant
Dimensional shrinkage	< 1%
Durability	Inert, stable
Reflectance	> 98%


Mechanical Design

(VII) Radiation Length

Radiation Length



Radiation length distribution of Aerogel Counter



✓ Overall radiation length is about 19.2% !!



- It is similar to TOF. (TOF ~ 18.7%)

Mechanical Design

(VIII) Aluminum Box




Box Production @ Dubna







Fabricated Box

✓ Aluminum Boxes were made in Dubna !!!

- Thickness of Lid : 0.8 mm
- Thickness of other parts : 0.5 mm



Assembled Counter




Apr 6th 2004

*PHENIX Focus ~ Aerogel Cherenkov Counter ~
Satoshi Takagi*

Electronics

(I) Front End Electronics

- ✓ ACC-FEE is the same as the RICH one, except for the trigger module.
 - 1 controller board, 5 AMU/ADC board
 - 1 readout board
- ✓ Read out signal from 160 PMTs
 - 0 to 50 p.e. detection with 10 bit
- ✓ AMU/ADC module
 - 64 channels for Analog to Digital Converter.
 - Charge and Timing information is digitized on AMU/ADC module.



Electronics

(II) Preamplifier

➤ Preamplifier Gain



- PMT signal : 0.8 pC/p.e.
- FEE dynamic range : 0 to 160 pC

To measure the photo-electrons at the range from 0 to 50 p.e. (0 to 40 pC)

➡ Required Preamp's net gain : $\times 4$

➤ Performance




- The measured gain is 4.3.
(design value is 4.5.)



Performance of Single Cell

Test Beam (KEK-PS)


(I) Clean Signal



- ✓ Very clean separation !!
- ✓ Amount of photons other than Areogel Cherenkov is small !!

Performance of Single Cell



Test Beam (KEK-PS) (II) Uniform Response




- ✓ Uniform response, thanks to Integration Volume
 - Important to separate in the momentum region of slow rise
- ✓ ~10% diff. between normal/reverse, due to diffusive nature of aerogel

Simulation Activity

(I) Optical Simulation



	<i>Sim.</i>	<i>Exp.</i>
<i>Np.e.</i>	<i>~ 14 p.e</i>	<i>~ 15p.e.</i>
<i>Propagation</i>	<i>~ 7 ns</i>	<i>~ 7 ns</i>
<i>Time diff.</i>	<i>~ 4ns (max)</i>	<i>~ 2 ns (max)</i>




Photon Propagation

Simulation Activity

(II) Aerogel Cherenkov Counter in PISA

Almost material put into PISA !!



Online Monitoring

(I) Information of the Summed PMTs


Charge distribution

- Convert ADC to number of p.e.
- The sum of p.e. in all PMTs. (160 PMTs).

Time distribution

- The time distribution of ACC in all PMTs(160 PMTs)


Aerogel Counter monitor
Run #: 110279 Date: Fri Jan 23 02:22:17 2004





Online monitoring

(II) Information of each PMT



Aerogel Counter monitor
Run #: 110279 Date: Fri Jan 23 02:22:17 2004



- Hit information
- Number of hits for each PMT



- Number of p.e.
- Number of mean p.e for each PMT



- Time information
- Mean TAC Value for each PMT

Calibration Method

➤ Conversion Parameter


- Conv. para. from ADC to p.e.
 - Our PMT has one p.e. resolution. (see right figure)
 - Fit the 1 p.e. peak by gaussian.

We calibrate ACC by using physics data.

➤ Pedestal Stability

- We took pedestal run once on a day.



We check pedestal stability by using daily pedestal run.



Zero-Suppressed data

Results @ Run4 Au+Au

(I) Raw data



- we can see the coincidence data clearly.
- Aerogel Counter works well !!

Results @ Run4 Au+Au


(II) Occupancy

➤ Simulation

- calculated by PISA
- HIJING Au+Au 200 GeV

➤ Experiment


- Run4 Au+Au 200 GeV




*Simulation result is consistent
with experimental result !!*

Results @ Run4 Au+Au

(III) Tracking association



Pion signal




Pion signal

- Clear pion rise up to $\sim 2\text{ GeV/c}$. (K, p not seen.)
- Peak position of Npe saturates at high momentum region.

Physics Results (future plan)

- Single particle p_T spectra
- Central-to-Peripheral ratio(R_{cp}) vs p_T
- Elliptic Flow
- etc...



In near future, the analysis results will be presented

Summary

- ✓ Half the detector has been installed in Run4.
 - remaining half detector will be installed by Run5.
 - New Time-of-Flight will be installed behind the Aerogel.
- ✓ The PHENIX Aerogel-Cherenkov-Counter is capable to extend the PID region of PHENIX.
 - pion identification : $\sim 3.7\text{GeV}/c$, $5.5\text{GeV}/c \sim 10\text{GeV}/c$
 - kaon identification : $\sim 3.7\text{GeV}/c$, $5.5\text{GeV}/c \sim 7\text{GeV}/c$
 - proton identification : $\sim 7\text{GeV}/c$
- ✓ In near future, the physics results will be presented.
 - example
 - Inclusive particle p_T spectra
 - Central-to-Peripheral ratio (R_{cp}) vs. p_T
 - Elliptic Flow
 - etc

Institutions

- ✓ BNL
- ✓ JINR-LHE (Dubna)
- ✓ CNS, University of Tokyo
- ✓ Tsukuba College of Technology
- ✓ University of Tsukuba

And, of course, much more people are related to the past/coming successful construction/operation of it.

Lots of persons for the essential works in the past, and much more persons in the future.