

Optical Alignment System for Muon Tracker

Atsushi Taketani,

RIKEN/RIKEN Brookhaven Research Center

Hiroki Kanou, Takashi Watanabe, Nobuyuki Kamihara,
and Takuma Horaguchi

Tokyo Institute of Technology

1. Hardware
2. Learned from RUN3/4
3. Upgrade
4. Plan and Summary

Optical Alignment System

Configuration

Light Source (station1)

- Single 150W Halogen light per
- Optical fiber to station1


Focusing lens(station2)

- 1cm convex lens

CCD camera (station3)

- effective area 8.8×6.6 mm
- Number of pixels 768×498
- Pixel size 11.0×13.0 μm


$$\text{Total 7optics/Octant} * 8 \text{ Octant/Arm} * 2\text{Arm} = 112$$



Data flow


CCD Camera

Getting spot image



DAQ PC

- Taking CCD camera image each 30 min.
- Making light intensity histogram for each X and Y projection.



Online PC

- Fit histogram and measure peak position



Position correction of Muon Tracking Chamber

- 1 . Measurement of peak position of light image.
- 2 . Make a model of motion of each Octant with parameters.
- 3 . Position correction
- 4 . Evaluation of correction.



Done by Takashi Watanabe as Master Thesis work

Measurement of peak position

Sharp image




Out focus Image



Long term movement

North octant4 CCD4



Long Term movement : 50~300 micron meter
on single CCD camera

Movement Model without expansion

Consider half octant as rigid body

$$\mathbf{f}_i = \begin{pmatrix} \mathbf{Rot}(\phi_x, \phi_y, \phi_z) \end{pmatrix} \begin{pmatrix} \mathbf{x}_i \\ \mathbf{y}_i \\ \mathbf{z}_i \end{pmatrix} + \begin{pmatrix} \Delta X \\ \Delta Y \\ \Delta Z \end{pmatrix}$$

$$\sim \begin{pmatrix} \text{Rotation} + \text{Position on Octant} + \text{Displacement in PHENIX} \end{pmatrix}$$

Look at the movement on X and Y

\mathbf{x}_i the position

Ignore rotation along X and Y axis

$$S = \text{Sum} |\mathbf{f}_i - \mathbf{x}_i|^2$$

Center of gravity of camera on octant

+ Δ

Peak position

+ Δ

$\Delta Y - \delta y_i)^2$

Minimize S with

in camera

placement

Movement model with expansion

$$\Delta \mathbf{r} = \alpha(t) \mathbf{r}: \text{Isotropic expansion}$$

$$\mathbf{f}_i = \begin{pmatrix} \text{Rot}(\phi_x, \phi_y, \phi_z) \end{pmatrix} (1 + \alpha) \begin{pmatrix} x_i \\ y_i \\ z_i \end{pmatrix} + \begin{pmatrix} \Delta X \\ \Delta Y \\ \Delta Z \end{pmatrix}$$

$$\mathbf{f}_i = \begin{pmatrix} (1 + \alpha)x_i - \phi_z y_i + \Delta X \\ \phi_z x_i + (1 + \alpha)y_i + \Delta Y \end{pmatrix} \quad : \text{Movement of Model}$$


$$\mathbf{x}_i = \begin{pmatrix} x_{0i} + \delta x_i \\ y_{0i} + \delta y_i \end{pmatrix} \quad : \text{Image position}$$

$$S = \text{Sum } |\mathbf{f}_i - \mathbf{x}_i|^2$$

Minimize S

$$= \text{Sum } \{ (\alpha x_{0i} - \phi_z y_{0i} + \Delta X - \delta x_i)^2 + (\phi_z x_{0i} + \alpha x_{0i} + \Delta Y - \delta y_i)^2 \}$$


Camera Combination



Combinatio	CCD4	CCD5	CCD6	CCD7
456	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
457	<input type="radio"/>	<input type="radio"/>		<input type="radio"/>
467	<input type="radio"/>		<input type="radio"/>	<input type="radio"/>
567		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

- Use North Arm, octant4, CCD4、 5、 6、 7
- Evaluate by comparing different camera configuration

Correction comparison on Δx



Correction evaluation

Displacement on X : 9~ 27 micron







Y : 5~ 24 micron

Rotation along Z axis : 4- 15 micro rad


Expansion : $2 \sim 15 \cdot 10^{-6}$

Temperature and Magnetic field dependence

- Temperature changes after magnetic field change
- Magnetic field moves chamber



RUN4



Data acquisition is not stable.

Pay attention by MuTr expert shift or PHENIX shift.

	North	South
Good	24	12
Recoverable	16	23
No Image	16	21

Good: good accuracy of peak

Recoverable : not good accuracy

No Image : No image at all

How to Improve

- No Image : Need to access the area inside magnet. -> Not this year.
- Recoverable: Take more picture, integrate them and then get sharper image.
- Replace DAQ system from GUI operation base to Labview base.

Daq system

Labview

Video Capture
board

Video signal

LAN

LAN

GPIB

Video
signal


Linux

Windows


GPIB-Ether Net

Multiplexer

CCD Camera
(*56)




Save data



Take pictures,
Integration
Make projection



Change cameras



All of components for prototyping are build at RIKEN

Executed 56 times {	Initialize system	few msec	
	Check hardware	few sec	
	Change camera	50 msec for 1ch	
	Take pictures	110 msec for 1ch, 1 picture	} Large contribution
	Integrate pictures	40 msec for 1ch, 1 picture	
	Make projection	40 msec for 1ch	
	Save to file	20 msec for 1ch	
	Close system	few msec	

Take N pictures for all camera...


$$T = 8.09N + 5.06 + \text{few (sec)}$$

($N < 73$, $T < 10$ min.)

Plan and Summary

- Young **Hiroki** will build Labview based system at RIKEN and move it to BNL at Aug.
- Looking OASYS will be shift duty.
- OASYS analysis data will be implemented in the geometry database and then improve J/PSI mass resolution.

Performance



Taking 1000 images from same camera for 30 minutes.

Temperature and Expansion

