

Performance of the PHENIX EMCalorimeter

QM2001 Poster

H.Torii, Kyoto Univ./RIKEN
for the PHENIX Collaboration

L. Aphecetche, T. Awes, S. Bathe, H. Baumsching, A. Bazilevsky, S. Belikov,
S. Botelho, S. Chernichenko, D. Bucher, H. Buesching, G. David,
A. Denisov, D. d'Enterria, H. Delagrang, A. Durum, R. Glasow, Y. Goto,
M. Grosse-Perdekamp, S. Fokin, H. Hamagaki, K. Imai, M. Ippolitov,
E. Kistenev, S. Kametani, V. Kochetkov, V. Manko, G. Martinez,
T. Matsumoto, Y. Melnikov, S. Mioduszewski, S. Nikolaev, A. Nyanin,
V. Onuchin, K. Oyama, Y. Schutz, A. Soldatov, T. Peitzmann, K. Reygers,
T. Sakaguchi, R. Santo, S. Stall, P. Stankus, T. Thomas, H. Torii, A. Usachev,
A. Vinogradov, M. Volkov, S. White, C. Woody, G. Young
for the PHENIX Collaboration

Physics Goal & EM probes



QGP search

- *Global characteristics*
- *Thermal radiation*
- *Chiral symmetry restoration*
- *Debye screening*
- *Jet quenching*
- *Space time evolution*
- *Fluctuation*

Spin physics

- ΔG
- $\Delta \bar{q}$

total transverse energy

direct γ , π^0

$\phi \rightarrow e^+e^-, K^+K^-$

$J/\Psi \rightarrow e^+e^-$

high- p_T π^0

HBT, flow

DCC

direct γ

π^0

open heavy flavor, J/Ψ

$W^\pm \rightarrow e^\pm$

PHENIX EM Calorimeter



- PbSc (lead scintillator) and PbGl (lead glass) -

Two type Calorimeter


- Lead Scintillator (PbSc)
 - 6 sectors (15552 channels)
- Lead Glass (PbGl)
 - 2 sectors (9216 channels)

Coverage

- $|\eta| < 0.38$
- $\phi = 180^\circ$

In the first year run


- $\phi = 45^\circ + 90^\circ$ are active.
 - 2 sectors PbSc
 - 1 sectors PbGl



PbSc

PbGl

PbSc Calorimeter



	PbSc
Size(cm x cm)	5.52 x 5.52
Depth(cm)	37.5
Number of towers	15552
Sampling fraction	~ 20%
h cov.	0.7
f cov.	90+45deg
h/mod	0.011
f/mod	0.011
X_0	18
Molière Radius	~ 3cm

Sandwich type calorimeter

Lead plates 55.2x55.2x1.5mm

Scintillator plates 110.4x110.4x4mm

Shish-kebab geometry wave shifter fiber readout

6x6 fibers $\rightarrow$ 1 PMT = 1 tower

2 x 2 towers = 1 module

6 x 6 module = 1 super module


6 x 3 super module = 1 sector



PbSc sector 2.0m x 4.0m

QM2001 poster

PbGl Calorimeter




	PbGl
Size(cm x cm)	4.0 x 4.0
Depth(cm)	40
Number of towers	9216
Sampling fraction	100%
h cov.	0.7
f cov.	45deg
h/mod	0.008
f/mod	0.008
X_0	14.4
Molière Radius	3.68cm

Lead Glass calorimeter

Lead Glass 40x40x400mm
used at WA98 exp.

4x6 towers = 1 super module

15*12 super module = 1 sector




PbGl sector 2.1m x 3.9m


Electron Measurement



E/p distribution for 0.3-1.0GeV/c electron, multiplicity<100




High multiplicity
condition



E/p distribution
shift up.

E/p multiplicity dependence




Electron E/p ratio is 1.0

The 'Ecorr' and 'Ecore' are the PHENIX EMCal clustering algorithm. 'Ecorr' is energy sum of all towers in a shower with linearity correction. 'Ecore' is energy sum of a few core towers to eliminate effect from adjacent shower overlapping. In high multiplicity condition, 'Ecore' shifts up by 3%, otherwise 'Ecorr' shifts by 7%.

Electron Energy measurement


~ Linearity ~



Test beam @AGS & CERN


Observed energy was corrected by the test beam result which is described by two different effects.

- Attenuation in fiber
- Shower Leakage



Electron Energy Measurement

~ Resolution ~




Test beam @AGS & CERN

Intrinsic energy resolution was obtained by test beam

$$1.9\% \oplus 8.2\%/\sqrt{E}$$

Current resolution is higher. It is under improvement.



Photon TOF measurement




Current status is

0.78nsec at $>0.3\text{GeV}$

0.67nsec at $>0.5\text{GeV}$

It is under improvement.


The intrinsic TOF resolution is measured at AGS beam test as 0.1nsec at $1\text{GeV}/c$.




π^0 measurement



Peripheral events




Subtraction by event mixing




75-92% centrality

Central events




Subtraction by event mixing




10% centrality

pT = 1.5 – 2.0 GeV/c Energy Asymmetry < 0.8

π^0 invariant mass




π^0 mass in PbSc



π^0 mass in PbGl


In higher multiplicity condition, the mass shows 1% higher than that in lower multiplicity condition. It is known to the overlapping effect of the two photon clusters. π^0 mass shows the energy calibration with accuracy of 2%.

Hadron PID




QM2001 poster

Hadron PID



0.6-0.8GeV/c



0.8-1.0GeV/c


π/K are separated up to 0.8GeV/c
Hadron physics will be obtained in
large acceptance by using low
energy hadron.

Hadron response

PHENIX

Energy deposit by charged hadron is characterized by.


- Range
 - π $\sim 0.3\text{GeV}/c$
 - K $\sim 0.5\text{GeV}/c$
 - p $\sim 0.9\text{GeV}/c$
- At larger momentum than range
 - Ionization energy
 - Caused by particles crossing the EMCal.
 - Minimum Ionization energy
 - Hadron shower.
- At less momentum than range
 - Annihilation process by anti-proton produce large energy observation.




MIP Energy



- # Comparison with
AGS test results: π 1.0GeV/c.
PHENIX $\sqrt{s_{NN}} = 130\text{GeV}$:
0.9-1.1GeV/c $\pi^\pm$ at lower multiplicity.
 - minimum Ionization Particle at 0.26GeV.
 - The spectrums of hadronic shower between 0.26 – 1.0GeV are very similar.
- # Relative energy of every towers are calibrated by using ionization energy of 0.6-1.0GeV/c $\pi^\pm$ at lower multiplicity. The accuracy of this calibration is 3% accuracy which is dominated by statistical error.



MIP Run dependence



Run dependence of minimum ionization energy is obtained from August through September. The stability is $< 2\%$

Summary



Electron measurement

- E/p distribution shows the peak at 1.0 from 0.3 to 1.0 GeV/c electron
- At higher multiplicity, E/p shifts by a few percent. It depends on the clustering algorithm.
- The width of E/p is wider than our expectation. It's under improvement.

γ/π^0 measurement

- The resolution of γ 's TOF is 780 psec at > 0.3 GeV. There is a room to be improved.
- 2-photon invariant mass spectrum shows π^0 peak clearly even at high multiplicity and high background condition.
- The measured π^0 mass is 137 MeV at peripheral events. Absolute energy calibration is consistent with 2% error.

Hadron measurement

- Hadron PID shows π/K separation up to 0.8 GeV/c.
- Observed energy is characterized for different particle and in different way over/under range.
- Relative energy was calibrated by using energy of minimum ionization particle at 0.26 GeV. The error is 3% dominated by statistical error.
- The comparison with the AGS beam test are very similar.
- The energy of minimum ionization particle is stable during the runs.